

Tilt Plate
QUASAR
BY MODUS QA



User's Guide

Version 2.0.0

MODUS QA
TRUSTED ACCURACY

Legal & Copyright Notice

Information in this manual is subject to change without notice. No part of this manual may be reproduced in any form without the written permission of Modus Medical Devices Inc (Modus QA).

Copyright ©2020 Modus Medical Devices Inc. All rights reserved. QUASAR™ and the QUASAR™ logo are trademarks of Modus Medical Devices Inc.

Modus QA reserves the right to make any changes without further notice to any products herein. Modus QA makes no representation, warranty or guarantee regarding the suitability of its products for any particular purpose, nor does Modus QA assume any liability arising out of the application or use of any product, and specifically disclaims any and all liability, including without limitation consequential or incidental damages. "Typical" parameters can and do vary in different applications. All operating parameters and default values, including "Typicals" must be validated for each customer application by customer's technical experts. Modus QA does not convey in this user's guide any license under its patent rights nor the rights of others.

Contents

Verify Image-Guided Shifts for 6° of Freedom Couches	5
Daily 6DoF QA.....	5
Specifications	5
Instructions.....	5
Daily 6DoF CBCT/Couch QA	5
Standard Daily CBCT/Couch QA	6
Additional Points	7
Appendix - Penta-Guide Tilt-Plate Dimensional References	8
<i>Anterior View</i>	8
<i>Lateral View</i>	8
<i>Axial View</i>	9
<i>Vertical Pillar</i>	9
Warranty.....	10
Modus Medical Devices Inc. Limited Warranty	10
Technical Support.....	10

This Page Intentionally Left Blank

Verify Image-Guided Shifts for 6° of Freedom Couches

Introducing QUASAR™ Penta-Guide Tilt Plate, an easy and efficient way to verify couch corrections for six degrees of freedom (6DoF) couches. The QUASAR™ Penta-Guide Tilt Plate positions the QUASAR™ Penta-Guide Phantom at known rotation angles and displacements relative to the isocenter. With this device, users test rotation corrections and combined translation and rotation corrections on a 6DoF couch.

Daily 6DoF QA

The central lines on the tilt plate, in each group of three, align the centre of the phantom with the laser isocenter for rotation only corrections. Other lines require translation and rotation corrections.

Specifications

- Tilt angles are: 0.75°, 1.25°, and 1.0° (see appendix for further dimensional references)
- Turning the tilt plate in 90° increments varies the direction of the rotation corrections
- Compatible with all new and existing QUASAR™ Penta-Guide Phantoms
- Compatible with most localization bars for faster, simpler test setup

Instructions

Daily QA of a linear accelerator includes testing of the treatment couch motion, as well as the on-board cone beam CT (CBCT) imaging system. Typically, this procedure is performed with the QUASAR™ Penta-Guide Phantom, which enables simultaneous testing of the couch motion (in 3 directions) and imaging systems. However, some treatment couches (eg. Hexapod Couch), are capable of 6 degree of freedom motion (3 translations and 3 rotations). Therefore, the QUASAR™ Penta-Guide Tilt Plate was designed as an accessory to the QUASAR™ Penta-Guide Phantom to facilitate the daily QA of linear accelerators equipped with these 6 degree of freedom couches.

Daily 6DoF CBCT/Couch QA

(Phantom with Tilt Plate for rotation only or rotation plus translation corrections)

1. Place the QUASAR™ Penta-Guide Phantom in the Tilt Plate
2. Align the Tilt Plate markings with the room lasers. The level indicators show that the Tilt Plate is level, but the phantom is not level.
3. Acquire CBCT of the QUASAR™ Penta-Guide Phantom and Tilt Plate

4. Use the image registration tools to align the CBCT with a reference CT of the QUASAR™ Penta-Guide Phantom (reference CT should not contain Tilt Plate) Ensure the calculated shifts are within 1mm of the expected shifts (magnitude of shifts depends on markings used), and that rotations are <math><0.5</math> degrees from expected rotation (- 5. Apply the correction and move the couch. Check that the lasers now coincide with the central markings on the QUASAR™ Penta-Guide Phantom. The level indicators show that the phantom is level, but the Tilt Plate is not level.




Figure 1 – Before correction: Lasers aligned with Tilt Plate




Figure 2 – After correction: Lasers aligned with Phantom

Standard Daily CBCT/Couch QA

(Phantom alone without the Tilt Plate for translation only corrections)

1. Align the off-centre markings on the QUASAR™ Penta-Guide Phantom with the lasers
2. Acquire CBCT of the QUASAR™ Penta-Guide Phantom
3. Use the image registration tools to align the CBCT with a reference CT of the QUASAR™ Penta-Guide Phantom (isocentre at the centre of the phantom on the reference CT)
4. Ensure the calculated shifts are within 1mm of the expected shifts (1cm, 1.2cm, and 1.4cm), and that rotations are <math><0.5</math> degrees (since there should be no rotations of the phantom)
5. Apply the correction and move the couch. Check that the lasers now coincide with the central markings on the QUASAR™ Penta-Guide Phantom

Additional Points

- Daily Couch/CBCT QA of a 6 degree of freedom couch with the tilt plate is simple. Only 1 extra step from the standard QUASAR™ Penta-Guide Phantom.
- When aligning the CBCT to the reference CT, it is best to perform an automatic match. However, the clipbox must be selected such that it does not contain the tilt plate. (Since there is no tilt plate in the reference scan).
- There are other vendor specific QC tests that need to be performed daily for the 6DoF couches - this does not replace those tests.
- There is a channel in the base of the Tilt Plate that is compatible with most localization bars
- Placing the Tilt Plate on a localization bar may simplify setup and improve repeatability of these tests

Appendix - Penta-Guide Tilt-Plate Dimensional References

Anterior View

Dimensions are in millimeters. Demonstrates that the central markings on the tilt-plate correspond to the center sphere location within the Penta-Guide phantom as it rests in the Tilt-Plate. The additional markings indicate a 10 mm displacement from that origin location. The yaw is indicated as a 1° counter-clockwise displacement.




Figure 3 – Anterior View

Lateral View

Indicating a 1.25° inclination (pitch) built into the Tilt-Plate surface, raising the superior aspect of the Penta-Guide anteriorly.




Figure 4 – Lateral View

Axial View

From the inferior end of the Tilt-Plate indicating a 0.75° inclination (roll) raising the Left (90° side) of the Penta-Guide phantom.




Figure 5 – Axial View

Vertical Pillar

The vertical pillar contains 3 reference lines. The central line aligns with the central sphere of the Penta-Guide phantom. The anterior and posterior lines indicate a 10mm displacement from the center line or center of the central Penta-Guide sphere.



Figure 6 – Vertical Pillar

Please note, these measurements are provided for reference only. As part of the recommended QA process, the residual error in alignment is the metric that should determine pass/fail criteria. I.e., after alignment images and corrections are applied, a verification image set (CBCT) should be performed to verify that the overall 6DoF alignment accuracy is within your expected tolerances (1-2mm).

Warranty

Modus Medical Devices Inc. Limited Warranty

Modus Medical Devices Inc. warrants to the original purchaser that the products sold will be free from defects in material and workmanship for a period of one year from the date of delivery. If the product is found to be defective during the warranty period, the part(s) and labor required to repair the product will be provided free of charge. This warranty is subject to the following exceptions and limitations:

- The customer shall be responsible for proper maintenance and handling of the product.
- No warranty is extended to any equipment that has been altered or modified in any way.
- No warranty is extended to any equipment that has been misused, or damaged.

This warranty is in lieu of all other warranties, express or implied, including warranties of merchantability or fitness for a particular purpose and is the sole and exclusive remedy for any claim of damages arising from defects in the Modus product. Modus shall have no liability for consequential damages or personal injury or for the loss, damage or expense, directly or indirectly arising from the use of this product. No distributor, dealer or other party is authorized to make any warranty on behalf of Modus, or to assume for Modus any other liability with respect to its products.

Technical Support

If you require service or support for this product, please contact Modus Medical Devices Inc. at:

Modus Medical Devices Inc.
1570 North Routledge Park
London, Ontario N6H 5L6
CANADA

Toll Free (North America): +1-866-862-9682
Tel: (519) 438-2409 x111

support@modusQA.com
www.modusQA.com



info@modusQA.com
www.modusQA.com



+1 (519) 438-2409
+1 (866) 862-9682 (NA)



Modus Medical Devices Inc.
1570 North Routledge Park
London, ON, Canada N6H 5L6

MODUS QA
TRUSTED ACCURACY